Your Original Tank-Mix Partner: ZeroTol® 2.0

⇔BioSafeSystems

Add to every spray for improved disease control

- Disease Reset Button
 On-contact activity resets the disease clock, reducing pathogen inoculum and pressure immediately
- **Enhances Efficacy of Systemic Chemistries**Allows focused activity of systemic or translaminar chemistries
- Preserves Fungicide Performance
 ZeroTol 2.0 helps prolong the life of standard chemistries and improves resistance management
- Cost-Effective Application Rates
 Rate ranges for every problem, every program, and every budget



Rates Explained

1:100

PRESSURE

High

Standalone applications only

1:300

PRESSURE

Med

Warmer months and sensitive crops

1:500

PRESSURE LOW

Proactive, low disease pressure months

Apply weekly at low rates, consistent use at 1:300 or 1:500 is more effective than intermittent and high-rate applications.

Do not mix with metal-based chemistries.